

## HIGH PRECISION ACCESSORIES HEAD SUPPORTS – LOW DENSITY

Models	Pictures	Dimensions
<p>Model 1 with lateral neck flaps</p> <p>AIO SOLUTION® code: <b>HS-NL-F</b></p> <p>Article No. : <b>35755</b></p>		 <p>100%</p>
<p>Model 2 without lateral neck flaps</p> <p>AIO SOLUTION® code: <b>HS-NL-N</b></p> <p>Article No. : <b>35765</b></p>		 <p>100%</p>
<p>Model 3 (age: 6-12 years) with lateral neck flaps</p> <p>AIO SOLUTION® code: <b>HS-PL-F</b></p> <p>Article No. : <b>35753</b></p>		 <p>90%</p>
<p>Model 4 (age: 6-12 years) without lateral neck flaps</p> <p>AIO SOLUTION® code: <b>HS-PL-N</b></p> <p>Article No. : <b>35753ZF</b></p>		 <p>90%</p>
<p>Model 5 with lateral neck flaps</p> <p>AIO SOLUTION® code: <b>HS-SL-F</b></p> <p>Article No. : <b>35752</b></p>		 <p>88%</p>
<p>Model 6 without lateral neck flaps</p> <p>AIO SOLUTION® code: <b>HS-SL-N</b></p> <p>Article No. : <b>35752ZF</b></p>		 <p>88%</p>
<p>Article No. : <b>32704</b></p>	<p>Set of Model 1, 2, 3, 4, 5 and 6</p>	
<p>Article No. : <b>32703</b></p>	<p>Set of 3 head supports with lateral neck flaps (Article Nos.: 35755 – 35753 – 35752)</p>	
<p>Article No. : <b>32705</b></p>	<p>Set of head supports without lateral neck flaps (Article Nos.: 35765 – 35753ZF – 35752ZF)</p>	

**A. GENERAL PRODUCT INFORMATION**

The products referred to in these instructions are medical devices, used for patient positioning and immobilisation in radiation therapy. These low density head supports can be used during both the simulation and the treatment stage, including MRI simulation.

To attain an optimal result, it is recommended to use this product in combination with Orfit immobilisation products.

**B. PRODUCT DESCRIPTION**

The HP HEAD SUPPORTS - LOW DENSITY are used in combination with the base plates, EFFICAST® PRE-CUTS, cushions and other hardware from the HIGH PRECISION and the AIO SOLUTION® sets. Together with the parts of these sets they form a reproducible patient positioning and immobilisation device in the field of radiotherapy. The head supports offer, together with the blocks & wedges, the opportunity to position the patient's head in a stable and reproducible way. Information on these other parts and instructions on how to make the masks can be found in the respective 'instructions for use' and on [www.orfit.com](http://www.orfit.com).

**C. PRECAUTIONS FOR USE**

Always use these devices in combination with HP or AIO base plates. The head supports can be positioned on the base plates or blocks & wedges by means of positioning knobs on the bottom of the head support. Always verify that the devices are positioned correctly on other cushions and parts.

This product is MRI-safe.



**D. STORAGE**

Always store the products in a safe place to prevent them from getting damaged. Do not put heavy objects on the devices and prevent hard objects from falling onto it, because these will cause permanent damage.

Store the devices between +10°C (50°F) and 40°C (122°F).

**E. PROPERTIES**

**E.1. Physical Properties**

The following physical properties apply to these devices:

Dimensions: L 230 - 257 mm  
W 108 - 120 mm  
H 82 - 92 mm  
Weight: 0.05 – 0.07 kg

**E.2. Mechanical Properties**

The LD head supports are made of PE foam covered with a PU coating resulting in a good dimensional stability. The cushions used on the treatment machines are often used more frequently than those on the simulators and during the imaging stages. We therefore recommend rotating the sets between the different machines on a regular base to avoid differences in wear and tear of the cushions.

The maximum deformation measured after 100 uninterrupted pressure cycles with a compression strength of 500N on a surface of 5 cm x 5 cm is 28.1%. This corresponds to the force of 100 times per day the load of 50 kg applied with the palm of the hand.

**E.3 Dosimetric Properties**

These head supports are made of low density foam that has the following dosimetric properties **per cm of material** through which the beam passes:

Attenuation (± 0.15 %)		Skin build-up (± 0.1 mm)
6 MV	15 MV	mm H <sub>2</sub> O equiv.
0.15 %	0.05 %	0.39 mm/cm material

Note: Use these numbers as a guidance only. Perform the measurements again in your department to verify these results.

**F. MAINTENANCE AND WASTE MANGEMENT**

These products can be cleaned and disinfected by means of an isopropanol based disinfectant, applied with a soft cloth. If unsure about the cleaning fluid, do not use. **Never use aerosol sprays, corrosive cleaning agents, solvents or abrasive detergents. Do not soak the cushions.** Further cleaning instructions can be found in the [Orfit Cleaning Guidelines](#). Periodic checks of the product should be done to insure the parts are not worn and require repair or replacement. **Do not attempt to make repairs yourself.** Contact your distributor if there are any questions or concerns. The products can be disposed of with household waste.

**G. ADDITIONAL INFORMATION**

For additional information such as distributor contact information, product brochures, Safety Data Sheets and regulatory information, please visit our website [www.orfit.com](http://www.orfit.com).

**Note:**  
It is prohibited to make alterations to this text without prior approval from Orfit Industries.  
THE AIO SOLUTION® and EFFICAST® are registered trademarks of Orfit Industries.



**ORFIT INDUSTRIES**  
Vosveld 9A | B-2110 Wijnegem | Belgium  
T (+32) (0)3 326 20 26  
welcome@orfit.com



Ref. No. 31826  
VERSION 20  
LAST UPDATE: 21/03/2023  
REVISION DATE: 21/03/2025